



mediland
Product Catalog

Company Profile

Established in 1984, Mediland is one of the major medical equipment manufacturers in Taiwan. It has been the Taiwanese and Chinese agent of GE Medical and has successful and rich sales experience in imaging medical equipment. Mediland was also the agent of worldwide leading medical equipment suppliers. With decades of professional medical equipment experience, Mediland has a complete distribution network in the medical industry.

Mediland entered the operating room equipment business in 1990 and manufactured products under the original brand Mediland. Major product ranges include surgical lights, surgical tables, surgical table accessories, operating room integration systems, disinfection system and pendant systems, etc. Mediland is a brand that fully represents the spirit of "made in Taiwan". At present, Mediland has production bases in Taoyuan, Taiwan and Nanjing, China.

Mediland brand value

Commitment

We are committed to taking the steps from simply being a supplier to becoming a valued partner and continuously improving the quality and technology of our products.

Service Orientation

As your trusted long-term partner, we listen to your requests and react quickly to tailor our service to your needs. We help our distributors to sell our products successfully in the global market resulting in profit for all parties.

Dedication

We are passionately dedicated to the development of complete OR solutions. Our experienced R&D team constantly improves the quality of our products and develops new high tech innovative solutions to meet our partner's needs. Our continuous research focuses on understanding changes in the market, bringing benefits to our clients and improving OR equipment for hospitals everywhere.



THE ONE YOU CAN TRUST



Surgical Table

Mediland's surgical tables deliver innovative technology and design to meet different surgical requirements. We offer different levels of tables with the same quality and unique features. For all surgical demands, Mediland is the brand you can rely on.

Amadis

Amadis operation table is built to meet the toughest conditions in today's OR to ensure the success of any surgical procedure. Amadis operation table allows Trendelenburg/Reverse Trendelenburg, lateral tilt and sliding adjustment.

- Hybrid power system for better stability
- High weight capacity suitable for different surgical procedure
- Wireless Hand Control for easy operation of the table
- Motorized drive wheel moves effortlessly in any direction



Artemis

Artemis is the table inherits unique design and reliable quality. This newly designed operation table portrays Mediland's innovative technology of high weight capacity. Artemis is the trustable partner in your high expectation operating environments.

- High weight capacity suitable for different surgical procedure
- Enlarged radiographic coverage distance
- Wireless Hand Control for easy operation of the table
- Aesthetic base design made of stainless steel resistant to cleaning substances used in the OR



Ares 8000

Ares 8000 is a high performance modular and motorized hydraulic surgical table independently developed by Mediland. The design incorporates innovative technologies to adapt to modern surgeries and allow for flexible patient positioning, thus improving the efficiency of medical procedures and ensuring the safety of patients.

- Innovative modular design
- High weight capacity
- Maximum Imaging Ares and sliding
- Optional Carbon Fiber construction
- Easy SmartLink System
- 400mm longitudinal Sliding



Amax 9000 & Amax9000 Plus

Amax 9000 Series is designed for better patient positioning and efficiency inside the OR. The table can adapt to all types of surgeries with functions such as Trendelenburg/Reverse Trendelenburg, lateral tilt kidney elevator and sliding adjustment.

- High weight capacity suitable for different surgical procedure
- Maximum range of coverage for C-Arm usage
- Table positioning can be powered by motor pump or foot pump
- T-shape base design provide sufficient legroom for surgeons



Surgical Table

SMART P2000

SMART P2000 is one unique model for its off-set column design for maximum X-ray coverage. 5-section tabletop design is flexible for various kinds of surgeries. Its slim column and compact base is perfect for sufficient legroom for surgeons.



- Rechargeable built-in batteries
- Off-set column design provides easy access for C-Arm
- 4 different models: Classic, Plus, Advantage and Elegance for choice to meet
- Different surgical requirement



C600

C600 is the most popular model for its universal application. Innovative design with multiple combination usage, it fulfills all operation position requirements.



- Rechargeable built-in batteries
- Optional features includes:
 - E: Extra shoulder section
 - K: Carbon fiber kidney elevator
 - S: Electric sliding up to 300mm
 - N: Low profile tabletop down to 520mm
 - D: Auto-leveling



C200

C200 has 3 different tabletop models, each configurable with up to 3 optional functions. These configurations provide various models for different markets and surgical preferences.



- Rechargeable built-in batteries
- Driven either electro-hydraulic or manual hydraulic system
- Optional functions includes :
S: Sliding
N: Low profile tabletop
B: Built-in batteries
- 300mm longitudinal sliding (S option)
- A robust and reliable control interface of the



Emax300

Emax300 electric surgical table can be used in the full range of surgeries. It is made of stainless steel capable of resisting all kinds of chemicals and disinfectants used in the OR. With its elegant design and reliable performance, the Emax300 will enhance the efficiency of your OR.



- Rechargeable built-in batteries
- Detachable leg section powered by gas springs
- Optional carbon fiber kidney elevator, leveling and sliding functions are available
- 4 different models: Classic, Plus, Advantage and Elegance for your choice to meet different surgical requirement



Surgical Light

MediLED Symphony

MediLED is the surgical lighting solution pioneering in the latest LED technology. It provides optimal visual environment while rendering accurate colors and appropriate light intensity. It delivers astonishing brilliance, longer lasting LED life span, more energy efficient and cooler temperature compared to traditional halogen solutions.

- Beautiful 7-petal formation around the central ring is good for laminar air flow compliant
- Easy to operate lightheads to desired position
- Brilliant nature white light with true colors
- Superb ability to render the accurate color of red by R9 index reaches 95
- Better shadow control
- User friendly touch screen, digital panel for light and camera control
- Three different lighthead modules to adjust to your OR requirements.
- Able for selection lightheads with intensities of 100,000, 140,000 and 160,000 Lux.
- Optional adjustable color temperature offers better enhancement and contrast on meticulous surgical site



SlimLED & SlimLED Plus

SlimLED series represents new lighting system from our LED surgical lighting product line. Surpassing traditional lighting systems in all aspects, this product brings perfect form, brilliant lighting, full functionality and cost savings to the ever demanding hospital administrations and surgeons.



- 6 cm compact, oval-shape design for more workspace for surgeons
- Easy to operate lightheads to desired position
- Brilliant nature white light with true colors
- Superb ability to render the accurate color of red by R9 index reaches 95
- LEDs distribution creates excellent shadow management control
- User friendly panel for light and camera control



ORLED

ORLED is the new generation of LED surgical light. It inherits the excellent performance and also combines variety of unique features and advantages. ORLED stands out as a high quality solution that will guarantee strong performance in any operation room.




- Patented double reflector design provides consistent homogeneous light column while also enhancing shadow control
- Extended flexibility in light field applications
- Easy to operate lightheads to desired position
- Special cooling system for heat dissipating
- Brilliant nature white light with true colors
- Superb ability to render the accurate color of red by R9 index reaches 97
- User friendly control panel for light
- Wireless in-light camera allows 360 degrees rotation for the lightheads




HyperLight Disinfection Robot

HyperLight Disinfection Robot Disinfection to Protection

Hyper Light Disinfection Robot is a specially designed robot made to prevent Healthcare-Associated Infection (HAI) for any healthcare facilities by means of UVGI disinfection method. The germicidal UV (254nm UVC) is capable of killing and inactivating microorganisms by destroying nucleic acids and disrupting DNA or RNA, leaving them unable to replicate and repair.



- Amalgam UV lamp features
- Rotational protective reflector technology
- Easy validation on Fast Check
- Intuitive wireless tablet control
- Multi-motion sensor system
- 3.2kg force conforming to AORN standards



ORPort

ORPort

ORPort Pendant System is used to integrate equipment and solve the space issue of overcrowding. The Pendant has a flexible moving range, equipped with various gas and power sockets, the system can provide the medical staff members a convenient and effective medical care environment.

The system is designed to be installed on the ceiling, this assists the pendant to be able to provide full rotation by the mechanical arms allowing the system to reach any desired location. ORPort Pendant System has a special exterior design, built from special aluminum alloy, its structure is safe and durable. ORPort Pendant System is configured with many accessories in order to meet the demands of users from different departments in the hospital. The firm is designed to bear large-size medical instruments making the ORPort an ideal medical gas, power supply, and instrument platform work station for medical staff members.





Mediland Enterprise Corporation

Linkou Factory

No.75, Wenhwa 1st Rd., Guishan Dist.,
Taoyuan City 33382, Taiwan, R.O.C.

Phone +886-3-328-9140
Fax +886-3-328-8881

international@mediland.com.tw
www.mediland.com.tw

Branch

No.29, Qinhuai Rd., Jiangning Economic &
Technological Development Zone, Nanjing,
Jiangsu, China

Phone +86-25-5791-8595
Fax +86-25-5791-8590



Printed in Taiwan, Jun. 2019 All rights reserved.
Mediland Catalogue Ref. No. QB-ProductCatalog-E-V2