

THE ONE YOU CAN TRUST



Surgical Table

Mediland's surgical tables deliver innovative technology and design to meet different surgical requirements. We offer different levels of tables with the same quality and unique features. For all surgical demands, Mediland is the brand you can rely on.

Amadis

Amadis operation table is built to meet the toughest conditions in today's OR to ensure the success of any surgical procedure. Amadis operation table allows Trendelenburg/Reverse Trendelenburg, lateral tilt and sliding adjustment.

- Hybrid power system for better stability
- High weight capacity suitable for different surgical procedure
- Wireless Hand Control for easy operation of the table
- Motorized drive wheel moves effortlessly in any direction



Artemis

Artemis is the table inherits unique design and reliable quality. This newly designed operation table portrays Mediland's innovative technology of high weight capacity. Artemis is the trustable partner in your high expectation operating environments.

- High weight capacity suitable for different surgical procedure
- Enlarged radiographic coverage distance
- Wireless Hand Control for easy operation of the table
- Aesthetic base design made of stainless steel resistant to cleaning substances used in the OR



Ares 8000

Ares 8000 is a high performance modular and motorized hydraulic surgical table independently developed by Mediland. The design incorporates innovative technologies to adapt to modern surgeries and allow for flexible patient positioning, thus improving the efficiency of medical procedures and ensuring the safety of patients.

- Innovative modular design
- High weight capacity
- Maximum Imaging Ares and sliding
- Optional Carbon Fiber construction
- Easy SmartLink System
- 400mm longitudinal Sliding



Amax 9000 & Amax9000 Plus

Amax 9000 Series is designed for better patient positioning and efficiency inside the OR. The table can adapt to all types of surgeries with functions such as Trendelenburg/Reverse Trendelenburg, lateral tilt kidney elevator and sliding adjustment.

- High weight capacity suitable for different surgical procedure
- Maximum range of coverage for C-Arm usage
- Table positioning can be powered by motor pump or foot pump
- T-shape base design provide sufficient legroom for surgeons



Surgical Table

SMART P2000

SMART P2000 is one unique model for its off-set column design for maximum X-ray coverage. 5-section tabletop design is flexible for various kinds of surgeries. Its slim column and compact base is perfect for sufficient legroom for surgeons.

- Rechargeable built-in batteries
- Off-set column design provides easy access for C-Arm
- 4 different models: Classic, Plus, Advantage and Elegance for choice to meet
- Different surgical requirement



C600

C600 is the most popular model for its universal application. Innovative design with multiple combination usage, it fulfills all operation position requirements.

- Rechargeable built-in batteries
- Optional features includes:
 - E: Extra shoulder section
 - K: Carbon fiber kidney elevator
 - S: Electric sliding up to 300mm
 - N: Low profile tabletop down to 520mm
 - D: Auto-leveling



C200

C200 has 3 different tabletop models, each configurable with up to 3 optional functions. These configurations provide various models for different markets and surgical preferences.



- Rechargeable built-in batteries
- Driven either electro-hydraulic or manual hydraulic system
- Optional functions includes :
S: Sliding
N: Low profile tabletop
B: Built-in batteries
- 300mm longitudinal sliding (S option)
- A robust and reliable control interface of the



Emax300

Emax300 electric surgical table can be used in the full range of surgeries. It is made of stainless steel capable of resisting all kinds of chemicals and disinfectants used in the OR. With its elegant design and reliable performance, the Emax300 will enhance the efficiency of your OR.



- Rechargeable built-in batteries
- Detachable leg section powered by gas springs
- Optional carbon fiber kidney elevator, leveling and sliding functions are available
- 4 different models: Classic, Plus, Advantage and Elegance for your choice to meet different surgical requirement

